

PHYSIOLITE II



oceanus

Precision Pain Relief



REGENERATIVE + TECHNOLOGY

66
My patients appreciate that I integrated the PhysioLite into their home care plan.

Dr. A. Volta, Denver, CO

LCD SCREEN

CONVENIENT SIZE

5 KEY CONTROL

AN EFFECTIVE THERAPY
FROM THE PRIVACY OF YOUR OWN HOME

Physical Effects of Pulsed Ultrasound

Your PhysioLite II is a therapeutic ultrasonic wave device that generates deeply penetrating heat that travels deep into the tissue and induce vasodilation, which helps increase blood flow to the treated area and is found to help restore your vitality faster.

Key Benefits

- > Clinically Proven Results
- > Accelerated Healing
- > Drug & Surgery Free

“


Amazing results after the first use! Faster recovery, less muscle soreness and im not sure if it is a side effect but my skin even looks better!!!

Abbey B.,
Minneapolis, Mn

”



How It Works



Treatment Recommendations

PhysioLite II therapy for soft tissue pain should be repeated daily to provide oxygen, nutrients and growth factors necessary for restoration.

1 SET DEVICE FOR TREATMENT

Program the ultrasound device with the provided settings:
Power Level, Duty Cycle, and Time

Energy 0.5-8 Watts | Duty Cycle 10-100% | Timer 15 Minutes

2 APPLY GEL

Place ultrasound gel (or approved alternative) to treatment areas.

3 START THERAPY

Press the ultrasound metal applicator against your treatment area and slowly glide across entire region for the duration.

Live Life Vibrantly+

TECHNICAL SPECIFICATIONS

- > **SonicEnergy: 1.0MHz**
- > **Pulsed Frequency: 100Hz**
- > **Duty Cycle 10% to 100%**
- > **Waveform: Continuous/Pulsed**
- > **Effective Radial Area: 5cm²**
- > **Max. Effective Intensity 1.6W/cm²**
- > **Preheat Time: Less Than 2 min.**
- > **Treatment Duration: <30 min.**
- > **FDA Product Code: ISA**

Clinical References

Aiyer R, Noori SA, Chang KV, et al. Therapeutic Ultrasound for Chronic Pain Management in Joints: A Systematic Review [published online ahead of print, 2019 May 16]. Pain Med. 2019;pnz102. doi:10.1093/pm/pnz102

**oceanus**

Oceanus America | Sales & Support 727-222-8573 | USA
www.oceanusamerica.com

Individual results may vary. Neither Oceanus nor any of its subsidiaries dispense medical advice. The content of this catalog does not constitute FDA approval, medical, legal, or any other professional advice.